DK50 2V

for two or three dental units



Technical specification	DK50 2V
Nominal voltage (V/Hz)	230 / 50 230 / 60 115 / 60 3x400/50
Maximal current of compressor (A)	7.5 8.5 16 4.5
Compressor efficiency (l.min ⁻¹ /6 bar)	135
Air tank capacity (I)	25
Pressure range (bar)	6 - 8
Mode of operation of compressor	S1 - 100%
Noise level (dB[A])	≤71
Weight (kg)	52
Dimensions (w x l x h) (mm)	460x460x690
Technical modifications are subject to change.	













DK50 2V/M

for two or three dental units

DK50 2V

without soundproof cabinet •

DK50 2VS

with soundproof cabinet •

DK50 2V/M


without cabinet, with membrane dryer •

DK50 2VS/M

with soundproof cabinet • with membrane dryer



Technical specification	DK50 2V/M	DK50 2V S/M
Nominal voltage (V/Hz)	230 / 50 230 / 60 115 / 60 3x400/50	
Maximal current of compressor (A)	7.8 8.8 16.3 4.7	8 9 16.4 4.9
Compressor efficiency (l.min ⁻¹ /6 bar)	110*	
Air tank capacity (I)	25	
Pressure range (bar)	6 - 8	
Mode of operation of compressor	S1 - 100%	
Drying point of compressor	Pressure dew point to - 20°C	
Noise level (dB[A])	≤71	≤ 53
Weight (kg)	57	94
Dimensions (w x x h)(mm)	460x535x690	560x665x860















DK50 2V/M

for two or three dental units















Technical specification	DK50 2V/M
Nominal voltage (V/Hz)	230/50
Maximal current of compressor (A)	7.8
Compressor efficiency (I.min ⁻¹ /6 bar)	90*
Air tank capacity (I)	25
Pressure range (bar)	6 - 8
Mode of operation of compressor	S1 - 100%
Drying point of compressor	Pressure dew point to - 40°C
Noise level (dB[A])	≤71
Weight (kg)	66
Dimensions (w x l x h) (mm)	460x535x690

* Part of the efficiency of compressor is used for regeneration of the dryer. Technical modifications are subject to change.



Absolute Air and Gas Ltd. Unit 6, Isidore Road Bromsgrove Technology Park Bromsgrove, B60 3FJ Phone: +44 845 605 6660 Fax: +44 1527 559262 enquiries@absolute air and gas. co.ukwww.absoluteairandgas.co.uk